


# Operating Instructions

## METTLER TOLEDO DL36 KF Coulometer





<b>Contents</b>	<b>Page</b>
<b>1. Introduction .....</b>	<b>3</b>
<b>2. Safety measures .....</b>	<b>4</b>
<b>3. Measurement principle .....</b>	<b>5</b>
<b>4. Putting into operation .....</b>	<b>6</b>
4.1 Preparing the titration cell .....	6
4.2 Attaching peripheral devices .....	7
4.3 Switching on the instrument .....	7
<b>5. Keypad and display .....</b>	<b>8</b>
5.1 The keypad .....	8
5.2 The display .....	9
<b>6. The first titration .....</b>	<b>10</b>
<b>7. The method .....</b>	<b>11</b>
7.1 Standard method .....	11
7.2 Printing out a method .....	12
7.3 Modifying a method .....	12
7.3.1 Titration parameters .....	12
7.3.2 Result parameters .....	14
<b>8. Setup .....</b>	<b>17</b>
8.1 Checking the reagent capacity .....	17
8.2 Entering the date and time .....	18
8.3 Entering the user identification .....	19
8.4 Instrument identification .....	19
8.5 Selecting peripheral devices .....	20
8.6 Selecting weight entry .....	21
8.7 Selecting entry of the sample identification .....	22
8.8 Initializing default parameters .....	22
<b>9. Sample determinations .....</b>	<b>23</b>
9.1 Weight transfer from a balance .....	23
9.2 Tips for sample addition .....	24
9.3 Tips for accurate determinations .....	25
9.4 Applications .....	25

---

	Page
<b>10. Error messages and malfunctions .....</b>	<b>26</b>
10.1 Error messages of the DL36 .....	26
10.2 Other malfunctions .....	26
<b>11. Maintenance and servicing .....</b>	<b>28</b>
11.1 Charging internal battery .....	28
11.2 Line voltage, line fuses .....	28
11.3 Greasing the glass joints .....	29
11.4 Changing the septum .....	29
11.5 Changing the silica gel .....	29
11.6 Changing the reagents .....	29
11.7 Cleaning the detection electrode .....	30
11.8 Cleaning and drying the inner burette .....	30
11.9 Measures for long-term storage .....	31
<b>12. Standard and optional equipment .....</b>	<b>32</b>
12.1 Standard equipment .....	32
12.2 Optional equipment .....	33
<b>13. Technical data .....</b>	<b>36</b>
13.1 DL36 .....	36
13.2 Interfaces of the DL36 .....	37
<b>14. Index .....</b>	<b>39</b>

## 1. Introduction

Moisture in sugar or salt influences its pouring properties. Traces of moisture in brake fluid have an adverse effect on the functioning of vehicle brakes. Impurity water in the hydraulic oil of aircraft leads to a dangerous lowering in the performance of the hydraulic system. As the quality of products frequently depends on their moisture content, it is essential that this be known as accurately as possible. The DL36 KF Coulometer can be used to determine trace moisture with high accuracy.

The DL36

- can determine the moisture content of a sample in a range from 10 µg to 100 mg H<sub>2</sub>O,
- allows entry of the sample weight before or after the titration either manually or automatically via an attached balance,
- automatically determines the drift and takes it into account in the result calculation,
- permits a recalculation of the results.

An attached printer prints out defined parameters and records the titration data in short-form or GLP format. An attached computer can interchange data with the titrator.

## 2. Safety measures

The DL36 has been tested for the experiments and intended purposes documented in the Operating Instructions. However, this does not absolve you from the responsibility of performing your own tests of the product supplied by us regarding its suitability for the methods and purposes you intend to use it for. You should therefore observe the following safety measures.

### Measures for your protection



- Ensure that you plug the power cable supplied into a receptacle outlet that is grounded! In the absence of grounding, a technical fault could be lethal.
- Switch the instrument off and disconnect the power cable before you change blown fuses! An electric shock could be lethal.



- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).
- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

### Measures for operational safety



- Check the set operating voltage before you switch on the instrument! The instrument may suffer damage if the operating voltage does not match the line voltage.
- Use only fuses of the type specified in the Operating Instructions! (Risk of fire)
- Have the instrument serviced only by METTLER TOLEDO Service!
- Always wipe off splashed liquids immediately! The instrument is not water-proof.
- Exclude the following environmental influences:
  - powerful vibrations,
  - direct sunlight,
  - atmospheric humidity greater than 80%,
  - temperatures below 5 °C and above 35 °C,
  - powerful electric or magnetic fields!

### 3. Measurement principle


In the Karl Fischer method for determination of the moisture content, water ( $\text{H}_2\text{O}$ ) reacts with iodine ( $\text{I}_2$ ) and sulfur dioxide ( $\text{SO}_2$ ) in the presence of methanol ( $\text{CH}_3\text{OH}$ ) and a base (RN).



In the volumetric titration, iodine is added as the titrant. In the coulometric titration, iodine is generated electrolytically by an iodide-containing anolyte.



As long as water is present in the titration cell, the generated iodine reacts according to reaction (1). As soon as all the water has been consumed by the reaction, there is a small excess of iodine in the anolyte. The double-pin platinum electrode detects this iodine excess and the iodine generation is stopped. According to Faraday's law, the amount of iodine generated is proportional to the current which has flowed. In reaction (1),  $\text{I}_2$  and  $\text{H}_2\text{O}$  react with each other in proportion 1:1. One mole of water (18 g) thus corresponds to  $2 \times 96\,500$  coulomb, in other words per mg  $\text{H}_2\text{O}$  a quantity of electricity of 10.72 coulomb is consumed. The total current consumption is a measure of the amount of moisture present.



## 4. Putting into operation

### 4.1 Preparing the titration cell


#### WARNING

Never inhale the vapors of Karl Fischer reagents and avoid skin contact! The reagents could hazard your health.

- Remove the Teflon stopper of the titration cell and carefully slip the magnetic stirrer into the cell; replace the stopper.
- Place the titration cell in the holder.

When the DL36 is delivered, the glass joints for the double-pin platinum electrode, inner burette and Teflon stopper have a very thin coating of silicon grease to assure tightness.


- Remove the Teflon stopper and add anolyte up to the lower mark on the titration cell (100 mL) using the funnel.
- Remove the silicon stopper from the inner burette and add 5 mL catholyte to the cathode of the inner burette; its level should be just **below** that of the anolyte.
- Check that the glass joint is greased, then insert the drying tube.
- Plug the cables of the electrode and the inner burette into the appropriate inputs at the rear of the DL36 (see illustration on next page).





## 4.2 Attaching peripheral devices

You can attach, e.g. the GA42 Printer and an analytical balance from METTLER TOLEDO and a computer (see Section 8.5: *Selecting peripheral devices* and Section 12: *Standard and optional equipment*).



## 4.3 Switching on the instrument

The DL36 operates in a voltage range of 100 - 120 or 220 - 240 V.



- Check the operating voltage set on the voltage selector before switching on the instrument! (See Section 11: *Maintenance and servicing*.)

- Plug in the power cable, connect to the power supply and switch on the instrument. The following appears in the display:

DL36 V1.0


followed by

< Run >

This display flashes until you press the **Run** key.

## 5. Keypad and display

### 5.1 The keypad



When you press this key, you can activate the stirrer or change the stirring speed; entry values: 0 to 9; with 0 the stirrer is at a standstill. The change is effected immediately.


- Press this key again to redisplay the current status.




If you press this key twice in quick succession, the stored method and setup parameters are printed out. You can print out the titration, result and setup parameters separately if you press the appropriate key and then press **Print** (see Sections 7.2 and 8).



When you stop the pretitration, <Run> appears in the display; when you stop the titration, e.g. "112 \* Pre-Titr" appears (see next page).



If you press this key after an entry/selection before you confirm the parameter with **Enter**, it will not be stored.



You must press this key if you

- wish to change the default parameters and values of the method/setup. The selection possibility is marked by an arrow: ↑ (see Section 7.1).
- wish to transfer the weight from an attached balance before or after the titration (see Section 9.1).
- wish to view the potential of the electrode during the pretitration and titration; pressing this key again switches back to the previous display (toggle key).

Pot.	201.5 mV
------	----------

Example with stable drift, this also means that the electrode potential must lie in the vicinity of this value at the end point of a titration.

## 5.2 The display

The display has space for 16 characters. Depending on the key selection, menu items, parameters or the prompt <Run> for the start of the titration are displayed: The following are shown during the titration:

- the drift in  $\mu\text{g}/\text{min}$
- the status \* character during pretitration and the titration (sample determination)  
      $\Rightarrow$  character for the start of a titration with the information: "Ready" or "Stable" (see next section)
- the changing, calculated moisture content.

112 * Pre-Titr
----------------

current drift value  
during the pretitration

479 *	206.2 $\mu\text{g}$
-------	---------------------

current drift  
value

amount of water titrated  
up to this point  
during the titration

## 6. The first titration

To determine the moisture content of a sample accurately, parameters are stored in the DL36 which you can change to match your sample. On delivery, the instrument has default parameters stored which you can use to determine the moisture content in  $\mu\text{g}$  of many organic substances, e.g. 99% ethanol.

The DL36 is switched on, the display shows **<Run>**.

– Press the **Run** key:

112 \* Pre-Titr

The pretitration starts; the solution is stirred and titrated until it is free from moisture before you inject the sample.

5  $\Rightarrow$  Ready

An audio signal sounds, this means the solution is anhydrous (ready); however, the drift is not yet stable, its value is  $\leq 500 \mu\text{g}/\text{min}$ . The drift is the amount of water titrated as moisture per minute and recorded. For accurate determinations, you should wait until

2  $\Rightarrow$  Stable

the drift is stable: it now changes by  $\leq 0.1 \mu\text{g}$  per minute.

– Prepare the sample: siphon off approx. 0.5 mL 99% ethanol ( $\text{H}_2\text{O}$  content = 0.2%) with a syringe.

– Press the **Run** key:

ID ? ETHANOL

– Enter the sample identification using the "minus", "period" and numeric keys (max. 10 characters, see Section 8.3). Confirm the entry with **Enter**:

2  $\Rightarrow$  Start  
Samp.In.

"Start" and "Samp.In." (sample injection) appear alternately in the display.

– Inject the sample and press the **Run** key:

978 \* 467.2  $\mu\text{g}$

The DL36 starts the titration and shows the constantly changing drift and calculated amount of moisture.

3 786.6  $\mu\text{g}$

When the titration is at an end, an audio signal sounds. The drift and the result are displayed and the result printed out on an attached printer.

- If you wish to determine the second sample immediately, prepare this and press the **Run** key: the sample identification will be requested (see above).
- When you press the **Enter** key, "Stable" (or "Ready") is displayed: the moisture diffusing into the cell will be continuously titrated, in other words the drift value is constantly redetermined (so-called standby titration).
- When you press the **Stop** key, "<Run>" is displayed.

## 7. The method

A method with default parameters is stored in the DL36. These comprise two groups, the measurement and result parameters, which you call up with the **Titration** and **Result** keys and can adapt to meet the requirements of the individual determinations. You can call up the parameters which support the method using the **Setup** key (see Section 8).

### Notes

1. During the titration, the **Titration**, **Result** and **Setup** keys are locked!
2. The stirring speed is not stored as a method parameter and is also not recorded. The stirrer is automatically activated when you start the pretitration (default value = 4). However, you can also activate the stirrer beforehand using the **Stirrer** key. If you switch off the instrument, the stirring speed you last defined remains stored.

### 7.1 Standard method

Group	Menu	Parameter	Value
Titration		t(stir) [s] ? t(wait) [s] ? t(max) [s] ? Drift stop: rel ↑ rel [µg/min] ?	0 15 0 6
Result	Calculation ↑	Unit: µg ↑ Weight: fix ↑ Weight? Drift comp: off ↑ Blank? Factor	5.0000 0 1.00
	Report ↑	Short ↑	
	Recalculation ↑	Weight ? Drift ? Blank ? Unit µg ↑	5.0000 0.00 0

?: A numeric entry or confirmation of the displayed value is expected.

↑: Menus or parameters have been preset which you can select with the **Select** key. You must then confirm the selection with **Enter**.

## 7.2 Printing out a method

You can always print out the parameters of a group:

- When you have selected **Titration** and press the **Print** key, the corresponding parameters are printed out.
- When you have selected **Result** and press the **Print** key, the calculation and report parameters are printed out. These parameters are printed out separately when you press the **Print** key while a parameter of the corresponding menu is displayed.

## 7.3 Modifying a method

You can change the parameters of the groups before and during the pretitration or standby titration.

### 7.3.1 Titration parameters

The parameters are responsible for the control of the titration. You define stirring and titration times as well as drift values which lead to abort of the titration.

Key selection	Display/Entry (example)	Description
<input type="text" value="Titration"/>	t(stir) [s]? 0 ..... <b>10</b>	Stirring time after the start of the titration: used to ensure thorough mixing or complete release of moisture from the sample; during this time, <b>no iodine is generated</b> , in other words no H <sub>2</sub> O is titrated.
<input type="text" value="Enter"/>	t(wait) [s]? 15	Wait time: the titration can not be aborted during this time. A value <b>&lt;15</b> is not possible!
<input type="text" value="Enter"/>	t(max) [s]? 0	t(max) is an abort parameter, the titration is aborted after "t(max) + t(wait)" even if the end point has not been reached (t(max) = <b>0</b> : abort parameter not defined, see next page).
<input type="text" value="Enter"/>	Drift stop: rel ↑	The drift stop is the second abort parameter. rel: relative drift stop; the titration is aborted when the drift is less than the "measured drift (before the titration) + the defined value" of, e.g. 4 µg/min.
<input type="text" value="Select"/>	Drift stop: abs ↑	abs: absolute drift stop; the titration is aborted when the drift is less than the entered drift value, e.g. 6 µg/min.
<input type="text" value="Select"/>	Drift stop: off ↑	off: abort parameter not defined.
<input type="text" value="Select"/>	Drift stop: rel ↑	<b>Notice:</b> You must define either t(max) > <b>0</b> or a drift stop value, otherwise you receive an error message during the titration. If you define both parameters, the titration will be aborted as soon as one of the values is reached.
<input type="text" value="Enter"/>	rel [µg//min] 6 ..... <b>4</b>	Enter, e.g. <b>4</b> for the relative drift stop.
<input type="text" value="Enter"/>		"<Run>", "Ready" or "Stable" is displayed.

## Notes

### 1. Stirring and titration time

Run	Run		
Sample injection	Stirring time $t(\text{stir})$	$t(\text{wait})$ : no abort of the titration	$t(\text{max})$ or the time until defined drift stop is reached
No iodine generation		Iodine generation	
Automatic drift compensation by internal calculation		Titration time for which you can select the type of drift compensation (see page 14)	

The titration time is printed out (Titr. time) if you select "GLP" for the report (see Section 7.3.2).

### 2. Abort parameters $t(\text{max})$ and drift stop


A defined titration time " $t(\text{wait}) + t(\text{max})$ " is suitable for titrations in which the end point is reached only very slowly owing to side-reactions.

With  $t(\text{max}) = 0$ , the titration time depends on the drift stop.

A higher value for the drift stop shortens the duration of the titration, but leads to a larger titration error, above all when the sample releases moisture only slowly.

- With oil samples, for instance, an increase in the drift can be observed. In this case, the end point will not be reached if the drift stop value is too low.
- In the case of titrations with a sluggish approach to the end point as a result of side-reactions, the magnitude of the error can be reduced by a higher drift stop value.

Titration curves with rapid and slow release of moisture with two relative drift stop values:



A: measured drift before the titration  
 L1: high drift stop value  
 L2: low drift stop value

End point 1: with high drift stop value  
 End point 2: with low drift stop value

### 7.3.2 Result parameters

This parameter group comprises three parts:

1. the calculation parameters (Calculation)
2. the selection of the report (Report)
3. the parameters for the recalculation (Recalculation).

The tables shown follow this subdivision.

#### 1. Calculation parameters

Key selection	Display/Entry (example)	Description
<input type="button" value="Result"/>	Calculation ↑	Select the result unit for the calculation of the moisture content (the absolute µg content is always recorded). The choice of unit automatically determines the formula used for the calculation (see next page).
<input type="button" value="Enter"/>	Unit: µg ↑	
<input type="button" value="Select"/>	Unit: mg ↑	
<input type="button" value="Select"/>	Unit: ppm ↑	
<input type="button" value="Select"/>	Unit: % ↑	Select the entry of the weight [g]: var: variable weight; inquiry of the weight for every titration. fix: fixed value; the weight for the samples being titrated does not change, the inquiry does not appear again in the titration. If you need only the absolute moisture content without weight value, enter <b>0</b> (zero).
<input type="button" value="Enter"/>	Weight: var ↑	
<input type="button" value="Select"/>	Weight: fix ↑	
<input type="button" value="Enter"/>	Weight? 5.0000 ..... <b>0.1</b>	
<input type="button" value="Enter"/>	Drift comp: off ↑	Select the drift compensation for the titration: off: the drift value calculated before the titration is neither incorporated in the calculation nor recorded. manu: you enter a drift value which is incorporated in the calculation. auto: the drift value calculated before the titration is automatically incorporated in the calculation.
<input type="button" value="Select"/>	Drift comp: manu ↑	
<input type="button" value="Select"/>	Drift comp: auto ↑	
<input type="button" value="Enter"/>	Blank? 0 ..... <b>2.5</b>	Enter a blank value [µg], e.g. if you have dissolved a sample in a solvent whose moisture content you have determined beforehand. This blank value is taken into account in the calculation (see next page).
<input type="button" value="Enter"/>	Factor? 1.0	Enter a factor used to multiply the result, e.g. if you wish to convert the result to an amount of sample which you can not titrate.
<input type="button" value="Enter"/>	Calculation ↑	To titrate, press <b>Esc</b> , to select the report parameters, press <b>Select</b> .



## Calculation formulae

Formula No.	Result in	Formula
1	µg or mg	$\text{factor} \cdot (T_{\text{value}} - \text{drift} \cdot t - \text{blank})$
2	ppm or % (entry of the sample weight)	$\text{factor} \cdot \frac{(T_{\text{value}} - \text{drift} \cdot t - \text{blank})}{\text{weight}}$

$T_{\text{value}}$ : Result of the titration in µg

drift: Drift value in µg/min

t: Titration time in min (see page 13)

blank: Blank value in µg

weight: Sample weight in g

## 2. Report parameters

You can choose between two report formats: short-form and GLP format (GLP: Good Laboratory Practice).

Key selection	Display	Description
<input type="button" value="Result"/>	Calculation ↑	Select either the short-form or GLP format.
<input type="button" value="Select"/>	Report ↑	
<input type="button" value="Enter"/>	Short ↑	
<input type="button" value="Select"/>	GLP ↑	
<input type="button" value="Enter"/>	Report ↑	
		To titrate, press <b>Esc</b> .

### Short-form report

```
Date 96/06/02 13:26
Sample ID ETHANOL
Weight      0.1243 g
Result      219.9 µg
            0.17689 %
```

### GLP report

```
Date 96/06/02 13:26
Instr. ID DL36-82
Sample ID ETHANOL
Weight      0.1243 g
Factor       1.00
Result      219.9 µg
            0.17689 %
Titr.Time   00:00:26
Drift       1.62 µg/min
Blank       0 µg
User ID     KVE
```

### On printout:

```
— Date and time
— Instrument identification
— Sample identification
— Weight
— Factor
— Result
— Titration time
— Drift value
— Blank value
— User identification
```

### 3. Parameters for recalculation

If, e.g. you have made a mistake when entering the weight or have selected a different result unit, you can perform recalculations after every sample determination.

Key selection	Display/Entry (example)	Description
<input type="button" value="Result"/>	Calculation ↑	
<input type="button" value="Select"/>	Report ↑	
<input type="button" value="Select"/>	Recalculation ↑	
<input type="button" value="Enter"/>	Weight? 0.1243..... <b>0.1234</b>	Change the weight if need be.
<input type="button" value="Enter"/>	Drift? 1.62	Change the drift value if need be.
<input type="button" value="Enter"/>	Blank? 0	Enter a blank value if need be.
<input type="button" value="Enter"/>	Unit: µg ↑	Select a different result unit if need be, e.g. <b>ppm</b> .
<input type="button" value="Select"/>	Unit: mg ↑	
<input type="button" value="Select"/>	Unit: ppm ↑	
<input type="button" value="Select"/>	Unit: % ↑	
<input type="button" value="Enter"/>	1756.0	The result is displayed and printed out with the selected report parameters.

## 8. Setup

Under this key you call up menus and functions which support the method and allow selection of the peripheral devices.

When you press the key, the following appears

Setup[1-8]?

When you enter one of the numbers and confirm with **Enter**, you can define or select the following parameters:

- 1: Reagent capacity
- 2: Date & time
- 3: User identification (User ID)
- 4: Instrument identification (Instr ID)
- 5: Peripheral devices (Option)
- 6: Weight entry (Wt. entry)
- 7: Entry of the sample identification (ID entry)
- 8: Initialize default parameters (Memory clear).

Press the **Print** key if you wish to print out the stored data of the menus.

### 8.1 Checking the reagent capacity

The catholyte and anolyte solutions have a limited capacity for the moisture determination: The capacity limit of

- 100 mL anolyte is reached after the determination of approx. 1000 mg H<sub>2</sub>O.
- 5 mL catholyte is reached after the determination of approx. 300 mg H<sub>2</sub>O (see specifications of the reagents).

Exhausted reagents result in values of the moisture content which are too low and lead to longer titration times.

The DL36 calculates and sums the amount of moisture in mg which has been determined with anolyte and catholyte and displays these values. You can

- reset these values following reagent change,
- enter limit values for the capacity and
- have the DL36 inform you when these limit values are exceeded.

Key selection	Display/Entry (example)	Description
<b>Setup</b>	Setup [1-8]? <b>1</b>	"Reagent capacity" flashes on the display.
<b>Enter</b>	A.Capa. [mg]? 61	Display of the summed electrolytic current for the anolyte – converted to mg H <sub>2</sub> O. When you change the anolyte, enter <b>0</b> .
<b>Enter</b>	C.Capa. [mg]? 61	Display of the summed electrolytic current for the catholyte – converted to mg H <sub>2</sub> O. When you change the catholyte, enter <b>0</b> .
<b>Enter</b>	Alarm set: off ↑	You are not informed when the defined capacity limits are exceeded (see below).
<b>Select</b>	Alarm set: on ↑	You are informed when the defined capacity limits are exceeded (see below). After a titration, the display then shows "A.Capacity Over".
<b>Enter</b>	A.Alarm[mg]? 1000 .. <b>850</b>	If necessary, change the default value for the capacity limit of the anolyte.
<b>Enter</b>	C.Alarm[mg]? 300 ... <b>250</b>	If necessary, change the default value for the capacity limit of the catholyte.
<b>Enter</b>	Setup [1-8]?	Press <b>Esc</b> if you wish to have the current status displayed.

## 8.2 Entering the date and time

The DL36 has an internal clock. When the instrument is delivered, default values are set for the date and time (95/06/01 and 00:00) from which they will continuously update. You should change both values to ensure correct data are reported.

Key selection	Display/Entry (example)	Description
<b>Setup</b>	Setup [1-8]? <b>2</b>	"Date ? Time" flashes on the display.
<b>Enter</b>	YY MM DD? 950728 ..... <b>960718</b>	Enter year, month and day.
<b>Enter</b>	HH MM? 1051 ..... <b>1346</b>	Enter hours and minutes.
<b>Enter</b>	Setup [1-8]?	Press <b>Esc</b> if you wish to have the current status displayed.

### 8.3 Entering the user identification

In this menu you can enter your name; it appears on every GLP report.

You can enter the user and the instrument identification alphanumerically using the "minus", "period" and numeric keys.



When you press this key for the first time, the letter **A** is displayed, for the second time **B** → **C** → ... **Z** → . (period) → – (minus) → "space" → **A** →...



When you press this key for the first time, a "space" is displayed, for the second time a – (minus) → . (period) → **Z** → **Y** → ... **A** → – "space" →...

You must confirm each letter (character) with **Enter** for it to be accepted.

Key selection	Display/Entry (example)	Description
Setup	Setup [1-8]? 3	"User ID" flashes in the display.
Enter	A	Enter, e.g. "A. Aichert": Press the "-" key.
Enter	A.	Press the "." key repeatedly until a period appears.
Enter	A._	Press the "-" key repeatedly until a space appears etc.
Enter	A. _	You can enter maximum 15 characters.
•		
•		

### 8.4 Instrument identification

On delivery of the DL36, its correct identification is stored (see model plate at the rear); it appears on every GLP report. You can change the identification.

Key selection	Display/Entry (example)	Description
Setup	Setup [1-8]? 4	"Instr ID" flashes in the display.
Enter	MMA06408 ..... D	Enter, e.g. "DL36-1A": Press the "-" key repeatedly until <b>D</b> appears.
Enter	DL	Press the "-" key repeatedly until <b>L</b> appears.
Enter	DL3	Press the "3" key etc.
•		You can enter maximum 8 characters.
•		

## 8.5 Selecting peripheral devices

In this menu you select the printer and the balance you have attached to the DL36. For a computer you need to define the corresponding configuration data. The pin assignment of the RS232C interfaces is described in Section 13.2.

### 1. Printer

You can attach a METTLER TOLEDO GA42 or another commercially available printer. In the latter case, you need to define the configuration data.

Note: If you attach the GA42, you must set its DIP switch 2 to position ON (see GA42 Operating Instructions).

Key selection	Display/Entry	Description
<b>Setup</b>	Setup [1-8]? <b>5</b>	"Option" flashes on the display.
<b>Enter</b>	Printer ↑	If "RS232C" or "Balance" is displayed, use <b>Select</b> to select "Printer".
<b>Enter</b>	Printer: GA- ↑	Choose between the GA42 and another printer. If you have attached the GA42, confirm with <b>Enter</b> (its configuration is preset, see Section 13.2, point 1): "Setup..." is again displayed.
<b>Select</b>	Printer: Other ↑	
<b>Enter</b>	Baudrate: 1200 ↑	Select the baud rate: 300, 600, 1200, 2400, 4800 or 9600.
<b>Enter</b>	Parity: none ↑	Select the parity: even, odd or none.
<b>Enter</b>	Data bits: 8 ↑	Select the data bits: 7 or 8.
<b>Enter</b>	Stop bits: 2 ↑	Select the stop bits: 1 or 2.
<b>Enter</b>	Setup [1-8]?	

### 2. Balance

You can attach a balance from METTLER TOLEDO, A&D, Shimadzu or Sartorius. The configuration is preset for every balance when you select it (see Section 13.2).

Key selection	Display/Entry	Description
<b>Setup</b>	Setup [1-8]? <b>5</b>	"Option" flashes in the display.
<b>Enter</b>	Balance ↑	If "RS232C" or "Printer" is displayed, use <b>Select</b> to select "Balance".
<b>Enter</b>	Bal: Mettler ↑	Select the name of the company whose balance you have attached.
<b>Enter</b>	Setup [1-8]?	

## Notes

1. For all balances you must select the default setting with unidirectional transmission mode "Send Cont.".
2. If you attach a balance using the LC-RS9 cable, you must set positions 7, 3 and 4 (switches left, middle and right).

### 3. Computer

If you attach a computer, you must select the configuration of the RS232C interface.

Key selection	Display/Entry	Description
<b>Setup</b>	Setup [1-8]? <b>5</b>	"Option" flashes on the display.
<b>Enter</b>	RS232C ↑	If "Printer" or "Balance" is displayed, use <b>Select</b> to select "RS232C".
<b>Enter</b>	Baudrate: 4800 ↑	Select the baud rate: 9600, 300, 600, 1200, 2400 or 4800.
<b>Enter</b>	Parity: none ↑	Select the parity: even, odd or none.
<b>Enter</b>	Data bits: 8 ↑	Select the data bits: 7 or 8.
<b>Enter</b>	Stop bits: 1 ↑	Select the stop bits: 2 or 1.
<b>Enter</b>	Soft HS: off ↑	Select software handshake: on or off.
<b>Enter</b>	Setup [1-8]?	

You can use the installed computer to transfer data in both directions in a specific format. A detailed description of the communication between the DL36 and a computer can be found in a separate set of instructions.

## 8.6 Selecting weight entry

You can choose whether you enter the weight before or after the titration.

Key selection	Display/Entry	Description
<b>Setup</b>	Setup [1-8]? <b>6</b>	"Wt.entry" flashes in the display.
<b>Enter</b>	Wt.entry:before ↑	You enter the weight before the titration.
<b>Select</b>	Wt.entry:after ↑	You enter the weight after the titration.
<b>Enter</b>	Setup [1-8]?	

## 8.7 Selecting entry of the sample identification

You can choose whether you enter the sample identification before or after the titration or not at all.

Key selection	Display/Entry	Description
<div>Setup</div> <div>Enter</div> <div>Enter</div>	Setup [1-8]? <b>7</b>  ID entry:before ↑  Setup [1-8]?	"ID entry" flashes on the display.  Select between before, after and off. off: no inquiry of the sample identification during the titration, but the identification last entered will be automatically recorded.

## 8.8 Initializing default parameters

You can delete all entries you have made for the method: The default parameters will then be reentered by the DL36 (see Section 7.1). These include

- the stirring speed (default value = 4)
- the display of the amount of moisture in mg titrated up to this point (default value for catholyte and anolyte = 0)  
the defined limit values for anolyte and catholyte (A = 1000 mg, C = 300 mg)  
the message indicating these limit values have been exceeded does not appear: "Alarm set:off" (see Section 8.1)
- date and time (default values = 95/06/01 and 00:00, see Section 8.2)
- the entry of the weight and the sample number (default setting: entry **before** the titration, see Sections 8.6 and 8.7).

Use the following procedure to set the default parameters:

Key selection	Display/Entry	Description
<div>Setup</div> <div>Enter</div> <div>Select</div> <div>Enter</div>	Setup [1-8]? <b>8</b>  Memory clr. :off ↑ Memory clr. :on ↑ Setup [1-8]?	"Memory clr." flashes in the display.  The initialization is not possible.  The initialization is possible:  The confirmation is accompanied by an audio signal. <b>Notice:</b> Immediately switch the instrument off, wait five seconds then switch it on again! If you do not do this, undefined values will be entered for all parameters!

After the instrument has been switched on, "SRAM Init." is displayed: the default parameters are initialized.

– Press any key: "DL36 V1.0" appears followed by "<Run>".



## 9. Sample determinations

### 9.1 Weight transfer from a balance

You have selected a variable weight as a calculation parameter (see Section 7.3.2) and "before the titration" as the entry for the weight and sample identification (see Sections 8.6 and 8.7).

**Notice:** The mass unit of the balance must be set to **g**!

The analyte solution is anhydrous, its drift stable (see Section 6).

- Siphon off your sample with a syringe, place the syringe on the balance and tare to **0**.
- Press the **Run** key:

ID ?

- Enter the sample identification and confirm with **Enter**.

2 ⇒ Start  
Samp.In.

"Start" and "Samp.In." (sample injection) appear alternately in the display.

- Inject the sample and replace the empty syringe on the balance.
- Press the **Run** key:

Weight? 0.0

The weight is requested.

- Press the **Select** key:

Weight? 0.1873

The value shown on the balance is transferred.

- Press the **Run** key: The titration starts. An audio signal sounds when the titration is at an end; the result is displayed and printed out.

If you have selected "after the titration" as the entry for weight and sample identification, both parameters are requested when the titration is at an end. An audio signal sounds and the following is displayed:

Weight? 0.0

The weight is requested.

- Press the **Select** key:

Weight? 0.1873

The weight shown on the balance is transferred.

- Confirm with **Enter**:

ID ?

- Enter the sample identification and confirm with **Enter**.

2 0.02343 %

The result is displayed and printed out after the calculation.

## 9.2 Tips for sample addition

The titration cell of the DL36 is set up for the injection of samples.

- During sample preparation and sample addition, work as quickly as possible to keep the effect of atmospheric moisture to a minimum. As most samples are stored and transported under standard conditions of temperature and pressure, you should also weigh in and inject under these conditions.
  - You can cool readily volatile samples before sampling (notice: moisture could then condense on the syringe!), viscose samples can be warmed.
  - With hygroscopic liquids, you should never take samples from the surface as the moisture content is highest at this point.
  - Samples in which moisture is present in the form of a dispersion can possibly be treated in an ultrasonic bath before addition.
- Never remove the Teflon stopper to add solid samples! This would result in an excessive amount of moisture entering the cell.
  - For determination of the moisture content in solids, you can dissolve these in an appropriate solvent or extract them and then inject the solutions. You must determine the moisture content of the corresponding solvent/extraction agent beforehand and enter it as a blank value (see Section 7.3.2).
  - With solids for which no suitable solvent is available for the titration, e.g. polymers, moisture can be driven off only by heating. METTLER TOLEDO offers drying ovens for such situations which lead the resulting water vapor through tubing to the titration cell of the DL36, where it can be determined.
- Adjust the amount of sample to the expected moisture content to achieve short titration times and allow as many titrations as possible to be performed with the same reagent solutions. The following table provides an overview of the amount to be weighed in.

Moisture content		Sample weight	
50 -	100% .....	<0,01 g	
10 -	50% .....	0,02 →	0,01 g
1 -	10% .....	0,05 →	0,01 g
0,1 -	1% .....	0,1 →	0,01 g
0,01 -	0,1% .....	1,0 →	0,1 g
0,001 -	0,01% .....	5 →	1,0 g
0,0001 -	0,001% .....	10 →	5 g

not suitable for coulometry

### 9.3 Tips for accurate determinations

1. We advise always using the DL36 in the operating mode "Pretitration" (standby titration).
2. Always wait for a stable drift before you inject the sample. This allows determination of an exact value for the drift change per minute, which is incorporated in the calculation of the result.
3. Titrate the cell overnight when you determine samples whose absolute moisture content is less than 50  $\mu\text{g}$  (standby titration).
4. With a low moisture content, rinse the syringe with the sample before you perform the first titration.
5. Check the DL36 by measurements with a reference substance, e.g. Hydranal<sup>®</sup> check solution from Riedel-de Haën.
6. Change the reagents when
  - the anolyte solution reaches a volume of 150 mL (second mark on the titration cell),
  - the capacity of the catholyte or anolyte is exhausted (see Section 11.6),
  - the drift is too high ( $>10 \mu\text{g}/\text{min}$ ).
7. Always titrate a sample series with the same stirring speed.

### 9.4 Applications

We have developed applications for volumetric and coulometric titrations which can be used to determine moisture in solvents, petroleum products, plastics, foods and cosmetics. The applications include

- instructions for sample preparation
- advice for substances which cause side reactions in the titration
- literature references.

They are listed in Section 12.2.

## 10. Error messages and malfunctions

### 10.1 Error messages of the DL36

Error message	Cause	Measures
A.Capacity Over!	The summed electrolytic current (converted to mg H <sub>2</sub> O) of the anolyte is greater than its capacity	<ul style="list-style-type: none"> <li>• Change anolyte</li> <li>• Set capacity of the anolyte to <b>0</b> (see Section 8.1)</li> </ul>
C.Capacity Over!	The summed electrolytic current (converted to mg H <sub>2</sub> O) of the catholyte is greater than its capacity	<ul style="list-style-type: none"> <li>• Change catholyte</li> <li>• Set capacity of the catholyte to <b>0</b> (see Section 8.1)</li> </ul>
Current error!	There is no current flowing between anode and cathode	<ul style="list-style-type: none"> <li>• Check anolyte and catholyte</li> <li>• Check supply leads of cathode and anode</li> </ul>
Elect open!	The detection electrode is not working properly (break)	Plug in electrode correctly
Elect short!	The detection electrode is not working properly (short circuit)	Rectify short circuit (align platinum pins correctly)
Meas. Over!	Measurement range exceeded. The amount of moisture in this sample is greater than 100 mg H <sub>2</sub> O	Change catholyte, inject smaller amount of sample
Over Titr.!	Overtitrated: the anolyte contains an excess of iodine	<ul style="list-style-type: none"> <li>• Clean anode</li> <li>• Protect instrument against direct sunlight</li> </ul>
Para set miss!	Neither of the two abort parameters t(max) and drift stop is defined	Define either t(max), the drift stop or both parameters
Pre Amp Err-XX	Error in the preamplifier	Contact METTLER TOLEDO Service

### 10.2 Other malfunctions

Fault	Possible causes	Measures
Instrument can not be switched on	No line voltage Power cable not connected Fuses faulty	Check power cable, on/off switch, voltage selector and fuses
Malfunctions of the stirrer	Stirrer speed set to <b>0</b>	Check stirring speed (stirrer key and size of the magnetic stirrer (35 mm))

Fault	Possible causes	Measures
Drift too high	<p>Leaks in the system</p> <p>Capacity limits of the anolyte/catholyte reached or "wrong" reagent used</p> <p>Membrane wetted with moisture or cracked</p> <p>Anode touching the membrane</p>	<p>Change septum and drying agent, (apply very thin coating of grease to glass joints)</p> <p>Change reagents</p> <p>Clean inner burette or replace</p> <p>Insert anode adjustment element between membrane and anode and set a distance between 1.5 mm and 2 mm</p>
Reagent consumption too high	<p>Detection electrode contaminated or faulty</p> <p>Reagents exposed to direct sunlight</p> <p>Detection electrode contaminated or faulty</p> <p>Capacity limits of the anolyte/catholyte reached</p> <p>Deposits on the anode surface</p> <p>Anode touching the membrane</p>	<p>Clean/replace detection electrode</p> <p>Protect instrument against direct sunlight</p> <p>Clean/replace detection electrode</p> <p>Change reagents</p> <p>Clean inner burette</p> <p>Insert anode adjustment element between membrane and anode and set a distance between 1.5 mm and 2 mm</p>
End point not reached or reached too late	<p>Drift stop set to "off"</p> <p>Drift stop value too low</p> <p>Side reactions</p> <p>Capacity limits of the anolyte/catholyte reached</p> <p>Detection electrode contaminated or faulty</p>	<p>Set "drift stop" to "on" and enter value</p> <p>Enter a larger drift stop value</p> <p>Use a different KF reagent (e.g. for ketones)</p> <p>Change reagents</p> <p>Clean/replace detection electrode</p>
Poor reproducibility of the results	<p>Moisture content of the sample too high/low</p> <p>Drift too high</p> <p>Reagent consumption too high</p> <p>Drift stop value too high</p> <p>Capacity limits of the anolyte/catholyte reached</p>	<p>Adjust amount of sample (see Section 9.2)</p> <p>See "Drift too high".</p> <p>See "Reagent consumption too high"</p> <p>Enter lower drift stop value</p> <p>Change reagents</p>

## 11. Maintenance and servicing

### 11.1 Charging internal battery

A battery is responsible for the internal clock. Contact the METTLER TOLEDO Service when it is discharged, i.e. when the date and time are no longer correctly displayed.

### 11.2 Line voltage, line fuses


The DL36 can operate in different voltage ranges. The voltage selector is used to select the range.

Voltage range/supply	Setting on voltage selector	Fuses
100 V, 50/60 Hz	100 VAC	T3,15L250V
120 V, 50/60 Hz	120 VAC	
220 V, 50/60 Hz	220 VAC	T1,6L250V
230/240 V, 50/60 Hz	230/240 VAC	

#### Changing voltage/fuses



- Switch the instrument off and disconnect the power cable!
- Open the cover above the power cable socket using a flat screwdriver.



- Take out the drum of the voltage selector and reinsert horizontally with the desired voltage.
- Remove both fuse holders using a screwdriver and insert the appropriate fuses (see above).
- Reinsert the fuse holders (note direction of arrow) and replace the cover.

The same procedure applies to the replacement of blown fuses.

### 11.3 Greasing the glass joints

- Turn the stopper/electrode/drying tube of the titration cell about once a week to check their contact and to prevent them from sticking.
- Apply a very thin coating of grease if the contact (tightness) is no longer assured.

**Notice:** If an excessive amount of grease is applied, it could enter the titration cell and increase the drift owing to the moisture it contains!

### 11.4 Changing the septum

- Change the septum of the injection port at frequent intervals. As old septums can easily break, air can enter the titration cell and increase the drift.

### 11.5 Changing the silica gel

You should change the silica gel as soon as it starts turning pink.

- Unscrew the cover.
- Change the silica gel and screw the cover back on.

### 11.6 Changing the reagents

**WARNING**

Switch the DL36 off and place it in a fume hood! Never inhale the vapors of Karl Fischer reagents and avoid skin contact! The reagents could hazard your health.

#### 1. Catholyte

You should change the catholyte when

- the capacity limit is reached (see Section 8.1),
- the drift is too high ( $>10 \mu\text{g}/\text{min}$ ) or the status "Stable" is no longer reached,
- the membrane is contaminated.

Changing the catholyte at regular intervals avoids high drift values and measurement errors.

- Siphon off the catholyte using the polyethylene bottle.
- Add 5 mL catholyte using a syringe (up to lower mark on the titration cell).

**Notice:** After the change, select **Setup 1** and enter **0** for "C.Capacity". This sets the value of the summed electric current (converted to  $\text{mg H}_2\text{O}$ ) of the catholyte to zero (see Section 8.1).

## 2. Anolyte

You should change the anolyte when

- the capacity limit is reached (see Section 8.1),
  - the drift is too high ( $>10 \mu\text{g}/\text{min}$ ) or the status "Stable" is no longer reached,
  - its volume in the titration cell has reached 150 mL (upper marking on the titration cell); if the anolyte is diluted too highly, the sensitivity is lowered and the titrations take longer.
- Siphon off the anolyte using the polyethylene bottle.
  - Add 100 mL anolyte to the titration cell (up to the lower marking on the titration cell).

**Notice:** After the change, select **Setup 1** and enter **0** for "A.Capacity". This sets the value of the summed electric current (converted to mg  $\text{H}_2\text{O}$ ) of the anolyte to zero (see Section 8.1).

### 11.7 Cleaning the detection electrode

- If the platinum pins of the electrode are contaminated, first clean with concentrated nitric acid or carbon tetrachloride and then with methanol.

### 11.8 Cleaning and drying the inner burette

#### **WARNING**

Place the DL36 in a fume hood! Never inhale the vapors of Karl Fischer reagents and avoid skin contact! The reagents could damage your health.

#### 1. Cleaning with methanol

- Switch off the DL36.
- Remove the inner burette from the titration cell and allow the catholyte to flow out.
- Degrease the joint surfaces with methanol.
- Rinse out the cathode of the inner burette two or three times with methanol.
- Add approx. 10 mL methanol and place in a beaker.
- Fill the beaker with methanol until the methanol level in the cathode is reached.
- Allow to stand for around 30 minutes, then empty the inner burette and dry it (see point 3).



**2. Cleaning with chromic acid** (if foreign matter has collected on the membrane or the platinum surface).

- Perform the first 3 steps as described under point 1.
- Add approx. 10 mL chromic acid and place in a beaker.
- Fill the beaker with chromic acid until the level in the cathode is reached.
- Allow to stand overnight.
- Empty the inner burette and wash the cathode five or six times with pure water until the yellow color disappears (rinse chromic acid from the membrane).
- Rinse with methanol and dry.

Note: Chromic acid cleaning solution: dissolve approx 1.5 g potassium dichromate in 100 mL concentrated sulfuric acid.

**3. Drying**

- Place the inner burette in a vacuum drying oven for at least two hours (temperature no higher than 50 °C).
- If no vacuum drying oven is available, you can use a hair dryer (blow warm air over the membrane of the inner burette until it is dry) or a desiccator.

**11.9 Measures for long-term storage**

We recommend cleaning the titration cell, disassembling it, drying the parts and then storing them in a desiccator.

## 12. Standard and optional equipment

You can reorder every part with an order number; the numbers in brackets refer to the quantity supplied with the DL36.

### 12.1 Standard equipment

#### Order No.

DL36 KF Coulometer

1 set Operating Instructions

in accordance with your order

1 memo card

in accordance with your order

1 power cable


in accordance with your order

1 set fuses (3)

in accordance with your order

1 jar silicon grease

—




Titration cell, complete

51107410


**comprising**  
Silicon stopper

—




Titration cell

51107411




Inner burette

51107412




Double-pin platinum electrode

51107413




Teflon stopper with septum  
(injection port)

105017



Magnetic stirrer (1)  
length 35 mm

105053

**Order No.**


Drying tube (1)  
• filled with silica gel

105016




Suction bottle made of polyethylene (1)

105123




Plastic funnel (1)

105124



Septum (10), Ø: 12/2 mm      set of 10


105074



Anode adjustment element (1)

105011


## 12.2 Optional equipment



Fuses  
T3,15L250V      set of 3  
T1,6L250V      set of 3


54286

18560




Silicon grease

71300




Silica gel (1 kg)

105080

**Order No.**



Syringe 1 mL                      set of 100  
Syringe 10 mL                   set of 100

71492  
71482




Injection needle                set of 12  
80 mm x 1.2 mm

71483




METTLER TOLEDO printer

GA42




Cable for GA42 Printer

51190362




Cable for DPU-201GS printer  
from Seiko Instruments

51106171




Cable for InkJet printers from  
Canon, HP and EPSON  
(D-submin. 9-pin, female, D-submin. 25-pin, male)

51190363




Cable for AT balance

229029



LC-RS9 cable for AG, AB, PB,  
and PR balances

229065



Computer cable (RS232C)  
DTE  
(D-submin. 9-pin, female, D-submin. 9-pin, female)

51190362

**Order No.**

Operating Instructions	German	51709541
	English	51709542
	French	51709543
Memo card	German	51709585
	English	51709586
	French	51709587
METTLER TOLEDO brochures Karl Fischer applications for	Foods, beverages, cosmetics	724478
	Chemicals, solvents, petroleum products, plastics	724354
	Karl Fischer The method at a glance	724573

## 13. Technical data

### 13.1 DL36

Measurement method	coulometric titration following Karl Fischer
Control system	timing control by the microprocessor with constant current pulses
End point detection	by alternating current polarization
Status display for a titration	Ready: Titration possible for rough results Stable: Titration possible for exact results
End point display	by audio signal
Stirring	with magnetic stirrer, length 35 mm
Capacity of the titration cell	maximum 150 mL
Measurement range	10 µg to 100 mg H <sub>2</sub> O content of a sample
Resolution	0.1 µg
Repeatability	the relative standard deviation of the titration of 1 mg H <sub>2</sub> O in methanol is less than 0.3 %.
Measurement time	90 - 120 seconds for the titration of 1 mg H <sub>2</sub> O in methanol
Drift	automatic compensation for the result calculation selectable
Display of the amount of moisture	from 0.1 µg to 999999 µg (after the titration)
Entry range for	0 - 999 seconds
• stirring time t(stir)	15 - 999 seconds
• wait time t(wait)	0 - 9999 seconds
• t(max)	0 - 999 µg/min
• drift stop	
Calculation function	Calculation of the moisture content <ul style="list-style-type: none"> <li>• in µg or mg</li> <li>• in ppm or % with known weight</li> </ul>


Housing	Polyester
Titration cell	Pyrex®
Display	LCD, 16 characters
Dimensions	width: 260 mm, depth: 395 mm, height: 270 cm
Weight	approx. 6.4 kg
Power supply	100/120/220/230/240 VAC ±10% 50/60 Hz 50 VA
Fuses	T1,6L250V; T3,15L250V
Ambient conditions	<ul style="list-style-type: none"> <li>• use only indoors</li> <li>• height up to 2000 m</li> <li>• ambient temperature: +5 to +35 °C</li> <li>• maximum relative atmospheric humidity 80% for temperatures up to 31 °C, linearly decreasing down to 50% relative humidity at 40 °C</li> <li>• overvoltage category: II</li> <li>• pollution degree: 2</li> </ul>

## 13.2 Interfaces of the DL36

### 1. Printer interface

RS232C interface for the attachment of a GA42 Printer or for various commercial printers


Pin assignment



## 2. Balance interface

RS232C interface for the attachment of a balance with unidirectional transmission mode "Send Cont."

Pin assignment




Preset configuration on attachment of one of the possible balances from:

Selection Default	Mettler	Shimadzu	A & D	Sartorius
Baud rate	2400	1200	2400	1200
Parity	even	none	even	none
Data bits	7	8	7	8
Stop bits	1	1	1	1

## 3. Computer interface

RS232C interface for the attachment of a computer or terminal

Pin assignment





---

## 14. Index

? (question mark) 11

↑ (arrow) 11

**Abort parameters** 12

Accurate determinations (tips) 25

Anolyte 5, 6, 17, 30

Applications 25

**Balance** 7, 20

Battery (internal) 28

Blank 14

**Calculation formulae** 15

Calculation parameters 14

Capacity of anolyte 17

Capacity of catholyte 17

Catholyte 5, 6, 17, 29

Computer 7, 21

Connections 7

**Date** 18

Default parameters 11, 22

Detection electrode 5, 30

Display 9

Drift 10, 29, 30

Drift change 10, 25

Drift compensation 13, 14

Drift stop 12, 13

**Error messages** 26

Esc key 8

**Factor** 14

Faraday's law 5

Fuses (change) 28

**GA42 printer** 20

Glass joints 6, 29

GLP report 15

ID entry 22

Initialization (default parameters) 22

Inner burette 5, 30

Inputs 7

Interfaces 37

**Karl Fischer reagents** 6, 29

Keypad 8

**Long-term storage (titration cell)** 31

**Maintenance** 28

Malfunctions 26

Measurement principle 5

Memory clr. 22

Method 11

Minus key 19

**Option** 20

Optional equipment 33

**Period key** 19

Peripheral devices 7, 20

Potential (electrode) 8

Print key 8

Printer 7, 20

Putting into operation 6

**Ready** 10

Reagent capacity 17

Reagent change 29

Recalculation 16

Report parameters 15

Result key 11, 14

Result parameters 14

**Safety measures** 4

Sample addition (tips) 24

Sample determinations 23

Sample ID(entification) 15, 22

Select key 8  
Septum 29  
Servicing 28  
Setup 17  
Setup key 11  
Short-form report 15  
Silica gel 29  
Stable 10  
Standard equipment 32  
Standard method 11  
Stirrer key 8, 11  
Stirring speed 11  
Stirring time  $t(\text{stir})$  12, 13  
Stop key 8

$t(\text{max})$  12, 13  
 $t(\text{stir})$  12, 13  
 $t(\text{wait})$  12, 13  
Technical data 36  
Time 18  
Titration cell 5  
Titration key 11, 12  
Titration parameters 12  
Titration sequence 10  
Titrationszeit 13

**Units** 14  
User identification 19

**Voltage** 28  
Voltage selector 28

**Weight entry** 21  
Weight transfer 23



**To protect your METTLER TOLEDO product's future:  
METTLER TOLEDO Service assures the quality, measuring accuracy and  
preservation of value of all METTLER TOLEDO products for years to come.  
Please send for full details about our attractive terms of service.  
Thank you.**

Printed on 100% chlorine-free paper, for the sake of our environment.



\*P51709542\*

Subject to technical changes and to the availability  
of the accessories supplied with the instruments.

© Mettler-Toledo GmbH 1996 ME-51709542 Printed in Switzerland 9606/2.12

**Mettler-Toledo GmbH, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach, Tel. (01) 806 77 11, Fax (01) 806 73 50, Internet: <http://www.mt.com>**

**AT Mettler-Toledo Ges.m.b.H., A-1100 Wien**

Tel. (01) 604 19 80, Fax (01) 604 28 80

**BE n.v. Mettler-Toledo s.a., B-1651 Lot**

Tel. (02) 334 02 11, Fax (02) 378 16 65

**CN Mettler-Toledo (Shanghai) Ltd., Shanghai 200233**

Tel. (21) 6485 04 35, Fax (21) 6485 33 51

**DE Mettler-Toledo GmbH, D-35353 Giessen**

Tel. (0641) 50 70, Fax (0641) 507 128

**ES Mettler-Toledo S.A.E., E-08038 Barcelona**

Tel. (03) 223 22 22, Fax (03) 223 02 71

**HK Mettler-Toledo (HK) Ltd., Kowloon**

Tel. (02) 744 12 21, Fax (02) 744 68 78

**HU Mettler-Toledo, KFT, H-1173 Budapest**

Tel. (01) 257 98 89, Fax (01) 257 70 30

**JP Mettler-Toledo K.K., Yokohama 231**

Tel. (45) 633 53 50, Fax (45) 664 96 50

**MY Mettler-Toledo (M), Sdn. Bhd. 47301 Petaling Jaya**

Tel. (03) 703 27 73, Fax (03) 703 87 73

**NL Mettler-Toledo B.V., NL-4000 HA Tiel**

Tel. (0344) 638 363, Fax (0344) 638 390

**RU Mettler-Toledo C.I.S. AG, 10 1000 Moskau**

Tel. (95) 921 92 11, Fax (95) 921 63 53

**SG Mettler-Toledo (S) Pte. Ltd., Singapore 139944**

Tel. (07) 786 779, Fax (07) 764 904

**SK Mettler-Toledo, service, s.r.o., SK-83103 Bratislava**

Tel. (07) 525 21 70, Fax (07) 525 21 70

**TH Mettler-Toledo (Thailand), Bangkok 10320**

Tel. (02) 719 64 80, Fax (02) 719 64 79

**UK Mettler-Toledo Ltd., Leicester, LE4 1AW**

Tel. (0116) 235 70 70, Fax (0116) 236 63 99

**AU Mettler-Toledo Ltd., Port Melbourne, Victoria 3207**

Tel. (03) 9 646 45 51, Fax (03) 9 645 39 35

**CH Mettler-Toledo (Schweiz) AG, CH-8606 Greifensee**

Tel. (01) 944 45 45, Fax (01) 944 45 10

**CZ Mettler-Toledo, spol. s.r.o., CZ-12000 Praha 2**

Tel. (02) 25 15 55, Fax (02) 242 475 83

**DK Mettler-Toledo A/S, DK-2600 Glostrup**

Tel. (43) 270 800, Fax (43) 270 828

**FR Mettler-Toledo s.a., F-78222 Viroflay**

Tel. (01) 309 717 17, Fax (01) 309 716 16

**HR Mettler-Toledo, d.o.o., HR-10000 Zagreb**

Tel. (01) 230 41 47, Fax (01) 233 63 17

**IT Mettler-Toledo S.p.A., I-20026 Novate Milanese**

Tel. (02) 333 321, Fax (02) 356 29 73

**KR Mettler-Toledo (Korea) Ltd., Seoul (135-090)**

Tel. (02) 518 20 04, Fax (02) 518 08 13

**MY Mettler-Toledo (S.E.A.), 47301 Petaling Jaya**

Tel. (03) 704 17 73, Fax (03) 703 17 72

**PL Mettler-Toledo, Sp. z o.o., PL-02-929 Warszawa**

Tel. (22) 651 92 32, Fax (22) 651 71 72

**SE Mettler-Toledo AB, S-12008 Stockholm**

Tel. (08) 702 50 00, Fax (08) 642 45 62

**SG Mettler-Toledo (S.E.A.) Pte. Ltd., Singapore 139944**

Tel. (07) 786 779, Fax (07) 786 639

**SL Mettler-Toledo, d.o.o., SL-61111 Ljubljana**

Tel. (06) 112 35 764, Fax (06) 127 45 75

**TW Mettler-Toledo Pac Rim AG, Taipei**

Tel. (62) 579 59 55, Fax (62) 579 59 77

**US Mettler-Toledo, Inc., Hightstown, NJ 08520-0071**

Tel. (609) 448-3000, Fax (609) 586-5451